


Volatile Matter Crucible



Code	Capacity (cc)	Diameter (mm)	Depth (mm)	Crucible Estimated Weight (g)	Flanged Cover Estimated Weight (g)
GVM10	10mm	25mm	28mm	7	6
GVM15	15mm	31mm	31mm	11	8
GVM20	20mm	33mm	35mm	15	9
GVM25	25mm	36mm	35mm	19	10
GVM30	30mm	40mm	36mm	22	11
GVM40	40mm	42mm	42mm	30	13
GVM50	50mm	45mm	44mm	38	14