



Mould




Code	A	B	C	Estimated Weight (g)
CM1	30/31mm	6mm	39mm	18 or 30
CM2	31/33mm	6mm	41mm	19 or 32
CM3	32/34mm	6mm	41mm	21 or 36
CM4	32/38mm	6mm	45mm	22 or 36
CM5	35/37mm	6mm	44mm	25 or 48
CM6	38/41mm	6mm	49mm	29 or 48
CM7	40/42mm	6mm	49mm	34 or 48
CM8	37/40mm	6mm	49mm	28 or 48
CM13	32/36mm	6mm	45mm	40 or 50

Dimple Crucible




Code	A	B	C	D	E	Estimated Weight (g)
CC1	37mm	32mm	20mm	5mm	43mm	26
CC2	37mm	32mm	20mm	5mm	43mm	30

Flat Base Crucible



Code	A	B	C	D	E	Estimated Weight (g)
CC3	35mm	35mm	17mm	3.5mm	40mm	21
CC4	35mm	35mm	21mm	3.5mm	44mm	26
CC8	35mm	35mm	21mm	3mm	43mm	60
CC5-N	37mm	32mm	26mm		43mm	30
CC5WB	37mm	32mm	25mm		43mm	30

Flat Base Side Dimple Crucible



Code	A	B	C	D	E	Estimated Weight (g)
CC5	37mm	32mm	26mm	5mm	43mm	30