



xrFuse Crucible



| Code | A | B | C | D | E | Estimated Weight (g) |
|------|------|------|------|-----|------|----------------------|
| DC1 | 36mm | 38mm | 19mm | 6mm | 44mm | 30 |
| DC2 | 36mm | 38mm | 19mm | 6mm | 44mm | 30 |

xrFuse Mould



| Code | A | B | C | Estimated Weight (g) |
|------|---------|-----|------|----------------------|
| DM1 | 39/41mm | 4mm | 50mm | 60-100 |
| DM2 | 37/39mm | 4mm | 50mm | 60-100 |
| DM5 | 32/34mm | 4mm | 42mm | 60-100 |
| DM6 | 30/32mm | 4mm | 42mm | 60-100 |